


# Best<sup>TM</sup> ABT Molecular Imaging

TEAMBEST GLOBAL®



www.teambest.com



## **NEW!** Best Sub-Compact Model 200 Self-Shielded Cyclotron with Optional Second Chemistry Module (Patent Pending)

- Low energy compact system — can be placed next to PET/CT
- Integrated self-shielded cyclotron and chemistry module
- Easy to operate push-button graphic interface
- Capable of producing:  $^{18}\text{F}$ FDG,  $\text{Na}^{18}\text{F}$ ,  $^{18}\text{F}$ -MISO,  $^{18}\text{F}$ FLT,  $^{18}\text{F}$ -Choline,  $^{18}\text{F}$ -DOPA,  $^{18}\text{F}$ -PSMA,  $^{11}\text{C}$ ,  $^{13}\text{N}$ ,  $^{68}\text{Ga}$  and more!
- Ideal for Nuclear Cardiology and other Medical Applications



## **NEW!** Introducing a Novel Target for Best Model 200 Sub-Compact Self-Shielded Cyclotron (Patent Pending)

- Improved and easier clamp design
- Improved cooling
- Capable of double the beam current as original design
- Increased yield of Radiopharmaceuticals per run
- New material for better heat transfer

TeamBest Global Companies © 2022

Best ABT Molecular Imaging 3024 Topside Business Park Drive, Louisville, TN 37777 USA

Best Cyclotron Systems 7643 Fullerton Road, Springfield, Virginia 22153 USA

Best ABT tel: 865 982 0098 BCS tel: 703 451 2378 [www.bestabt.com](http://www.bestabt.com) [www.bestcyclotron.com](http://www.bestcyclotron.com)

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA