

## MOSFET Patient Dose Verification Systems for IMRT, IGRT/Tomotherapy, VMAT, TBI, Radiology & Brachytherapy

### Portable DOSIMETER

- Stand-alone dosimetry system
- Portable & compact with LCD display
- Optional data export via USB
- No software or computer necessary



### mobileMOSFET Wireless

- Dose obtained in Real-Time
- Simplifies dosimetry & minimizes QA time
- Wireless with Bluetooth technology
- Export dosimetry data to R&V system

### MOSFET Dosimeter Features:

- Isotropic ( $\pm 2\%$  for  $360^\circ$ )
- Active area of  $0.2 \times 0.2$  mm
- Pinpoint measurement for high gradient fields
- Same calibration factor for photons & electrons at MeV energies



TeamBest Companies © 2019

USA orders: CNMC Company tel: 615 391 3076 800 635 2662 sales@cnmcco.com

CAN orders: Best Medical Canada tel: 613 591 2100 877 668 6636 bmcinfo@teambest.com

[www.cnmcco.com](http://www.cnmcco.com) [www.bestmedicalcanada.com](http://www.bestmedicalcanada.com) [www.teambest.com](http://www.teambest.com)

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA