

# Best™ Cyclotron Systems

TEAMBEST GLOBAL



www.teambest.com


Best Cyclotron Systems offers a large family of cyclotron designs with a broad energy range from 1 MeV to 150 MeV that can be configured and optimized for different applications.



NEW Best Cyclotrons	1–3 MeV	Deuterons for materials analysis (Patent Pending)
	70–150 MeV	For Proton Therapy (Patent Pending)
	3–90 MeV	High current proton beams for neutron production and delivery (Patent Pending)
Best 15p Cyclotron	15 MeV	Proton only, capable of high current up to 1000 Micro Amps, for medical radioisotopes
Best 20u/25p Cyclotrons	20, 25–15 MeV	Proton only, capable of high current up to 1000 Micro Amps, for medical radioisotopes
Best 30u/35p Cyclotrons	30, 35–15 MeV	Proton only, capable of high current up to 1000 Micro Amps, for medical radioisotopes
Best 70p Cyclotron	70–35 MeV	Proton only, capable of high current up to 1000 Micro Amps, for medical radioisotopes
Best 150p Cyclotron	From 70 MeV up to 150 MeV (non-variable)	For all Medical Treatments including Benign and Malignant Tumors for Neurological, Eye, Head/Neck, Pediatric, Lung Cancers, Vascular/Cardiac/Stenosis /Ablation, etc. (Patent Pending)



Installation of Best™ 70 MeV Cyclotron at INFN, Legnaro, Italy



### COMING SOON!

Best Proton Therapy Cyclotron up to 150 MeV dedicated for proton therapy with two beam lines and two treatment rooms (Patent Pending)

TeamBest Companies © 2020

TeamBest ASIA Private Limited email: [asia@teambest.com](mailto:asia@teambest.com)

2/299A Plot No. 36, 5th Cross Street • V.G.P. Prabhu Nagar, Perumbakkam • Chennai – 600100 India

[www.teambestasia.in](http://www.teambestasia.in) [www.bestcyclotron.com](http://www.bestcyclotron.com) [www.teambest.com](http://www.teambest.com)

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA