

BestTM medical international

Best Theratronics GammaBeam Systems have been saving lives of cancer patients for over 60 years around the world!

TEAMBEST GLOBALTM



www.teambest.com




Best GammaBeam 500 Total Body Irradiator System


The Total Body Irradiator GammaBeam 500 is a teletherapy unit designed to produce a large fixed rectangular radiation field at an extended source-to-skin distance in order to deliver total body irradiation. The unit can also be used for research and dosimetry purposes. Features include: dose rate up to 40 cGy/min at 2.5 m in air, record & verify and imaging capabilities, graphical touch screen data entry interface and motorized vertical motion of the head.

BestTM Integrated Brachytherapy Solutions

Best Fiducial Markers



Best I-125 & Pd-103 Seeds



BestMick Applicator

BestMick 10/15/20/25 Cartridges



Best Prostate Stabilization Needle


Best Applicator Needle




Best HDR Remote Afterloader

Best Medical is the only company that makes custom seeds & strands to your exact specifications — shipped within 24 hours, 7 days a week, sterile & non-sterile!

Best Flexi Needles



Best Brachytherapy Kit



Best Syed Multi-Channel Vaginal Applicator

Best Medical International, Inc. 7643 Fullerton Road, Springfield, VA 22153 USA
tel: 703 451 2378 800 336 4970 www.besttotalsolutions.com www.teambest.com

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA

* Certain products shown are not available for sale currently.