

BestTM ABT Molecular Imaging

TeamBest Global Cyclotron installations list is available at:



<http://www.bestcyclotron.com/about.html>

Best BG-75 BioMarker Generator with Optional Second Chemistry Module (Patent Pending)

- Push button graphic interface
- Single or batch dose production
- Capable of producing: ^{18}F FDG, Na^{18}F , ^{18}F -MISO, ^{18}F LT, ^{18}F -Choline, ^{18}F -DOPA, ^{18}F -PSMA and more!
- Final dose delivery to syringe or vial
- Automated quality control testing
- Integrated self-shielded cyclotron and chemistry module
- Complete production lab in a 5 x 5 meter area

NEW! Best Sub-Compact Model 200 BioMarker Generator with Optional Second Chemistry Module (Patent Pending)

- Low Energy Compact System
- Self Shielding
- Easy to operate
- Very little maintenance
- Capable of producing the following in addition to BG-75 isotopes: ^{18}F , ^{11}C , ^{13}N and ^{68}Ga
- Ideal for Nuclear Cardiology and other Medical Applications





BestTM Cyclotron Systems

NEW! Best Compact Cyclotron 6–15 MeV Variable Energy (Patent Pending)

- 1–1000 μA extracted beam current
- Capable of producing the following isotopes: ^{18}F , ^{68}Ga , ^{89}Zr , $^{99\text{m}}\text{Tc}$, ^{11}C , ^{13}N , ^{15}O , ^{64}Cu , ^{67}Ga , ^{111}In , ^{124}I , ^{225}Ac , ^{103}Pd and more!
- Simultaneous ^{18}F and ^{68}Ga at high yield
- Up to 5×10^{13} neutrons per second from external target
- 21 stripping foils at each stripping port for 2 minute rapid change

TeamBest Global Companies © 2021



www.bestabt.com • www.bestcyclotron.com • www.bestproton.com

Best ABT tel: 865 982 0098 • BCS tel: 604 681 3327 • BPT tel: 703 451 2378

AFRICA | ASIA | EUROPE | LATIN AMERICA | MIDDLE EAST | NORTH AMERICA

TEAMBEST GLOBAL[®]



www.teambest.com