

Useful for:

GYN, Prostate, and other Brachytherapy (HDR/LDR) procedures in Urology and General Practice

Best® Universal Chair/Table

Model Number: UCT-1000

1. Removable Head Rest
2. Heavy Duty Moisture Proof Ergonomic Vinyl Cushions
3. Stainless Steel Bilateral Cushioned Knee Supports
4. Side Rails accept Standard US Type Clamps and are reversible to accept UK Type Clamps
5. Padded Variable Height Intravenous IV Arm Rest
6. Central Locking Brake
7. Cantilevered Design provides ample clearance below Seat Area
8. All Electronic Motorized Movements via Hand Held Controller (Pendant) or Foot Pedal
9. Heavy Duty Multi-directional Casters



Electronic

- Variable Height Control
- Variable Back Support
- Variable Foot Support
- Horizontal Tilt 12 Degrees Head or Foot

4-Function Hand-Held Pendant (Control)

Removable Stainless Steel IV Pole



Best® Universal Chair/Table: Model UCT-1000

Dimensions

Chair/Table in Sitting Position

Width (with Side Rails) 72.4 cm
Seat Height to Floor (MIN) 60.5 cm

Chair/Table in Horizontal Position

Width (with Side Rails) 72.4 cm
Length 212.6 cm
Min. Height to Floor 60.5 cm
Max. Height to Floor 95.5 cm
Max. Patient Weight 400 Lbs



Universal Urology Exam Chair/Table Model UCT-1000

- Heavy Duty Moisture Proof Ergonomic Vinyl Cushions
- Carbon Fiber Low Attenuation Seat
- (2) Multidirectional Cushioned Knee Supports
- (2) Safety Side Rails
- (4) Heavy Duty Multidirectional Hospital Quality Casters
- (2) Universal Accessory Rails (Adapts to all Stepper Stabilizer Systems)
- (1) Adjustable Intravenous Fluid IV Pole
- (1) Headrest Classic
- (6) Accessory Clamps
- Central Locking Brake System
- Electronic, Motor Driven, Variable Height, plus or minus 35 cm.
- Electronic Variable Back Support
- Electronic Variable Leg Support
- Electronic Variable Tilt Total 24 degrees Head to Foot
- (1) Hand Held Multiple Position Control Pendant
- (1) Variable Height Control Pedal
- (1) Battery Backup System
- Maximum Patient Weight, 400 Pounds.

