

Best[™] nomos[®]

VISIT US IN BOOTHS 3900 & 4000, SOUTH HALL ENTRANCE


Cyber Sonalis[®] Ultrasound System

Featuring SimulView[™] Ultrasound Technology



The **Best[™] Cyber Sonalis[®] Ultrasound Imaging System** provides superior visualization of HDR, LDR, RF or Cryo-surgical procedures. Our patented **SimulView[™] Technology** provides simultaneous “live” views of the prostate in both planes, thereby increasing treatment accuracy and precision.

“Live” views captured with SimulView[™] Technology using Cyber Sonalis[®] Ultrasound System:



Additional Features:

- ✓ Capable of supporting over 20 different probe configurations, this system can be used in almost every discipline
- ✓ 21.5" touch screen technology
- ✓ Superior HD image resolution for improved procedure accuracy, speed and physician confidence level
- ✓ Patented probe design
- ✓ Sagittal array provides for 140 mm length of view encompassing the bladder, prostate and perineal tissue
- ✓ Advanced modular software design provides for future upgrade path
- ✓ Convertible to desktop or portable without losing any functionality
- ✓ Advanced drawing and editing tools which include user-defined line widths and colors
- ✓ Dual battery backup built-in to protect from power loss
- ✓ Multiple brachytherapy grid options with user-definable grid symbols and colors



*Certain products shown are not available for sale currently.

Advanced keyboard technology, totally sealed and self-healing antibacterial keyboard with SensoFoil[™] Technology

21.5" touch screen with multi-point touch capabilities



Multiple options for Probe Holder positions and gel bottle cups

Electric lift column for a wide range of height options and rotating support column for easy positioning

Independent rotating/locking casters

Easily accessible ECG and computer I/O ports

Center mounted probe for maximum efficiency

Large rear handles for easy moving

Power cord input port easily accessible

Rear cable management hooks

Sonalis[®] Transducers

We offer probes for the most common types of procedures and imaging, all with the enhanced image quality and high resolution that you would expect from the **BEST** name in ultrasounds!

8L2A Linear Array	16L5 Linear Array	10EC4A Endocavity
12L5A Linear Array	8V3 Phased Array	10BP4 Bi-Plane
14L3 Linear Array	4V2A Phased Array	8BP4 Bi-Plane
15LW4 Linear Array	5C2A Curved Array	8TE3 Trans-esophageal
15LA Linear Array	9MC3 Curved Array	Pedoff
15L4A Linear Array	8EC4A Endocavity	16HL7 High Frequency Linear Array
	XY-BI-Plane Phased Array	

Best Nomos • One Best Drive • Pittsburgh, Pennsylvania 15202 USA
tel: 412 312 6700 800 70 NOMOS www.nomos.com www.teambest.com

TeamBest Global Companies ©2026. All rights reserved.